



50ml Pistol Syringe

Weiyou® 50ml pistol syringe is made of high quality aluminum casting body with chrome plated. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.



Product Description

Weiyou® Plastic-Steel Veterinary Syringes is high quality plastic (PC or TPX) with UV additive. Weiyou® Plastic-Steel Veterinary Syringes is mainly used animal veterinary anti-epidemic and treatment. OEM customization is acceptable. Vaccination to prevent disease for all kinds of animals. Weiyou® Plastic-Steel Veterinary Syringes is continuous injection and accurate scale. Free samples, and small quantity can be supplied.

Specification

Item	50ml pistol syringe
Specification	50ml
Scale	1-5ml
Weight	870g
Size	24.5*19*4cm
Parts	two gaskets, a glass pipe, a box syringe needle(12pcs/box)

Features

- Aluminum casting body with chrome plated
- Brass luer-lock with chrome plated
- With calibration
- Glass tube easy for cleaning
- Suitable for the inoculation of livestock
- Clean well and oil with silicone oil, regularly replace washer





www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com

















